

## Night Splints

### OrthoFast™ Night Splint

This Patented Night Splint is designed to provide the patient relief from chronic plantar faciitis and Achilles tendonitis, to treat foot drop, or to be used Post-surgically. The OrthoFast Night Splint provides static passive dorsiflexion and can be used with or without the provided wedge.

Item #	Size	Foot Length
DO3NI	Medium	8" - 10"
DO4NI	Large	10" - 12"

The OrthoFast Night Splint is:

- Lightweight
- Comfortable - Padded for comfort while sleeping
- Semi-Weightbearing - not a substitute for a walking boot
- Universal (no left or right)



### OrthoFast 2-in-1 Ankle Brace

One ankle product for Treatment • Rehabilitation • Prevention  
Stirrups + Figure 8 Straps

#### INDICATIONS

With Stirrups

Indicated for acute ankle sprain management: compression, immobilization and rehabilitation.

Without Stirrups

Can be used for chronic rehabilitation and prevention of ankle sprain as well as management of mild foot drop.



2-in-1 Sizes	Item #	Men's Sizes	Women's Sizes
SM	2N11	5 - 7	6 - 8
MD	2N12	8 - 9	9 - 10
LG	2N13	10 - 13	11 - 14
XL	2N14	14 - 15	15+

The Ortho Slipper is a night splint that functions by stretching the plantar fascia via the windlass effect as opposed to traditional night splints which work by maintaining the ankle in a dorsiflexed position. Traditional night splints tend to be cumbersome, expensive, uncomfortable and have a low compliance rate while the

### Ortho Slipper



DynaSlipper is comfortable with a high rate of patient compliance.

Most patients will find significant relief after the first night of use and many patients will achieve satisfaction in one week.

Stock#	Ortho Slipper Size Chart
NDXS	XS W: 5-6½
NDS	S W: 7-8, M: 5-7
NDM	M W: 8½ - 10, M: 7½ - 9
NDL	L W: 10½+, M: 9½- 11½
NDXL	XL M: 12+

### Foot Smart Arch Brace



This curved brace supports the plantar fascia reducing inflammation and relieving pain. With two adjustable Velcro straps, the arch brace fits comfortably under the arch and can be easily attached and secured in seconds. Wear against the skin or over hosiery or socks. Won't change shoe fit.

Style#:	Size(Shoe Size)
AB1033	A: M 7-8 1/2, W 6 - 10 1/2
AB1034	B: M 9 - 13

### Darco Medical Surgical Shoe

Utilizing its famous semi-rigid sole and built with a MetaShank forefoot support, this shoe provides rigid control and protection. A breathable mesh upper conforms to bony abnormalities of the foot and has a reinforced heel counter and padded heel collar to stabilize the rear foot and reduce slippage.



No. DMS, ea.

Men's: Sm., Med., Lg., XL.  
Women's: Sm., Med., Lg.

### Darco Ortho-Wedge Healing Shoes

Specifically designed to protect the forefoot by removing most of the pressure from the metatarsal head and digits. A natural for the treatment of diabetic ulcerations or as a post-op shoe after surgeries to correct hallux valgus, hammer toes or other osteotomies with or without pins.



No. DOW, ea.

Color: Black  
Sizes: Unisex - Sm., Med., Lg., XL.

### Darco A.P.B. Hi All Purpose Boot

Designed with a closed toe to provide all-season protection. A high top design and ankle strap completely eliminate heel slippage. Great after surgery and in the treatment of the diabetic and insensitive foot.



No. APB, ea.

Color: Black  
Sizes: Unisex - Sm., Med., Lg., XL.

### Vacuum Press & Vacuum Pump

This vacuum press features a unique perimeter air flow, oversized 17"x10" work area with 3" depth, built-in pressure relief valve, and more. The Gast Pump features both an air cleaner and an exhaust silencer.



Complete Unit	No. HURC, ea.
Vacuum Press Frame	No. HURF, ea.
Vacuum Pump	No. HURP, ea.
Replacement Bladder	No. HURB, ea.

### One Shaft Grinder

This Grinder features a 1/4 HP 2-speed motor with a 6" shaft height. Single shaft on right side with a replaceable dust bag for easy convenient cleaning on the left side. Chuck remover allows easy changing of sanding drums.



One Shaft Grinder  
Sanding Disks

No. A9130 ea.

40 grit (Coarse) 60 grit (Medium) 120 Grit (Fine) No. 9137 ea.